


Laser Induced Fluorescence Spectroscopy and Molecular Imaging as Tools for Tumor Detection *in vivo*

B. Ebert

Physikalisch-Techn. Bundesanstalt, Section 8.31, PTB, Berlin, Germany


Outline:

- Spectroscopic identification of malignant regions in the gastrointestinal tract
- Receptor-targeted fluorescence imaging of animals
- Demarcation of lymph nodes by fluorescence imaging
- Spatial resolution in fluorescence imaging
- Fluorescence reference material



Fluorescence imaging of lymph nodes

Experimental setup




Heme biosynthesis

Succinyl-CoA + Glycine



5 - Aminolevulinic acid (ALA)



2 mol ALA

Porphobilinogen


4 mol Porphobilinogen

Uroporphyrinogen III

Coproporphyrinogen

Protoporphyrinogen IX

Protoporphyrin IX
(PpIX)




M: CH_3
 V: $\text{CH}=\text{CH}_2$
 P: $\text{CH}_2-\text{CH}_2-\text{COOH}$

Charité

Fachbereich Medizinische Messtechnik
PTB
 2007


Fluorescence spectroscopy of tumors endoscopic view



Charité

Fachbereich Medizinische Meßtechnik
PTB
2007


Fluorescence spectroscopy of malignant tissue



$$R = I(635 \text{ nm}, 20 \text{ ns}) / I(599 \text{ nm}, 0 \text{ ns})$$



Fluorescence spectroscopy of lymph nodes

Comparison with histology (cumulative frequency)




Cancer Research, 2001, 61, 991-999

Experimental setup: Gated fluorescence imaging



Receptor targeted NIR- imaging of mouse xenografts with fluorescent ligands



nature biotechnology, 2001, • 19, 327 - 331

In vivo fluorescence imaging of tumors



SSTR2 Tumor, NIR38:
20nmol /kg body weight

Octreotid

R = dPheCys-Phe-dTrp-Lys-Thr-Cys- Thr



Optical molecular imaging of lymph nodes using a targeted vascular contrast agent



strong fluorescence in the liver consistent with the hepatobiliary elimination pathway

J Biomed Opt. 2005, 10(4):41205

18th June 2007

Workshop: High Brightness Laser Sources




Bayer Schering Pharma




Optical molecular imaging of lymph nodes using a targeted vascular contrast agent

$$\text{Contrast: } K = (I_{\text{lymph}} - I_{\text{muscle}}) / I_{\text{muscle}}$$





High sensitivity, requiring as little as 0.25 nmol dye per animal

MECA 79
Cyanine dye
Control
conjugate




Experimental setup: Gated fluorescence imaging






Depth resolved fluorescence imaging

Fluorescent rod at different depth in a scattering solution $\mu'_s = 14 \text{ cm}^{-1}$

depth: 1 mm




depth: 12 mm




Depth resolved fluorescence imaging

Fluorescent rod at different depth in a scattering solution $\mu'_s = 14 \text{ cm}^{-1}$





Influence of geometry



2.2 ns

2.85 ns

4.0 ns




Charité

Fachbereich Medizinische Messtechnik
PTB
2007


Fluorescence reference material and application

Fluorescence excitation spectra



Bayer Schering Pharma


Fluorescence intensity in dependence on concentration of cyanine dye molecules




Bayer Schering Pharma



Scattering phantom with NIR96007, 2% on glass spheres



Monitoring of inflammation of ankle joints using cyanine dyes as contrast agents



Bayer Schering Pharma


Charité



18th June 2007

Workshop: High Brightness Laser Sources

Increase of fluorescence intensity in the right ankle joint after i.v. application of a cyanine dye



J Fluoresc. 2005 15(3):337-62



Conclusion

- ★ **Small malignant regions (dysplasias) in the gastrointestinal tract can be identified.**
- ★ **Light sources with high pulse energy and low repetition rates are preferable to suppress ambient light**
- ★ **Successful targeting of endothelial surface-expressed molecules - leukocyte homing**
- ★ **Determination of depth of inclusions by time resolved fluorescence imaging**
- ★ **Fluorescent phantoms are needed, however, to match optical properties of biological tissue is difficult**